

# DENTAL UNIT

D V P R E M I E R E - M



# DV- PREMIERE - M

ACCURACY EFFICIENCY STABILITY  
RELIABILITY & USER FRIENDLY



Touch pad for user Friendly operations

Stable design for your ultimate comfort

Integrated design to ensure its trustworthiness

Innovative control system to ensure precision  
and effectiveness

### **Cuspidor Water Unit**

Switch between Induction and manual water collection  
Automatic pressure clean water system  
Automatic Bowl Rinsing

### **Assistant Element**

Equipped with safety protection for securer operation  
Securer and easier operation for your beneficial experience  
Assistant Element is easy to operate with trustworthy hands



### **OPERATION DENTAL LIGHT**

LED for shadowless operation light  
Color Temperature 4000-5000K ,Brightness 25000 Lux  
Stepless light adjustment Removable handle for easy disinfection  
Multi-axis rotation design to ensure that the light is irradiated at any Angle Longer life

### **WORKING STOOLS**

360° Swiveling backrest and seat  
Manual hand release lever  
Mounted on 5-castors



# **PATIENT DENTAL CHAIR**

Synchronized movement of chair offer boundless comfort for patient  
Further safety system under backrest and base plate ensures operation safety  
Chair safety system is present for dental chair and dental unit to work simultaneously  
30 degrees in horizontal rotation angle of dental chair provides more operational space for doctor . Right Armrest is present for efficient positioning of patient  
Programmable Chair positions & pre-set positions electric driven  
Backrest movement range 1-70 degrees  
Movement Range Head rest 150mm





## Dentist Element

Time and alarm Settings Priority function of first device Fault display  
 Air locking system for easy manipulation, Rinsing time, Linkage function  
 Personalized design with 9 memory display: Analog display of instrument speed  
 Three different user interfaces can be provided the tray can be enlarged as required  
 Innovative & Latest instrument arm with two sets ball bearing ensures extraordinary performance



### DENTAL HANDPIECE

Optical fiber brightness can be adjusted at any level from 1 to 16



### FOOT PEDAL

Lift and back motion of the chair  
 Free rotation of rotatable instruments  
 Operation light can be controlled via  
 Efficient foot pedal



### DISINFECTION SYSTEM

Automatic pipe disinfection system  
 Easier operation

# STANDARDS ACCESSORIES

## DENTIST ELEMENT

- 3-Way Syringe
- X-ray Film Viewer
- Small X-ray Film Viewer
- 3 outlets for turbine or air motor without Light

## STOOL

- Gas Lift Cylinder
- Five castors
- Height adjustable lever

## OPERATING LIGHT

- LED Light

## WATER UNIT

- Moveable Ceramic Cuspidor
- Water Bottle System
- Heat for cup Filler
- Venturi System

## DENTAL CHAIR

- Trendelenburg Movement
- Multidimensional Control
- Safety protection Switches
- Seamless Upholstery
- Headrest 2 articulated
- Left & Right Armse
- Pre-set Positions



## ASSISTANT ELEMENT

- Control Panel
- 3-Way Syringe
- Large Suction Outlet
- Small Suction Outlet
- Moveable Assistant Element

## OPTIONS

- Scaler
- Air Compressor
- Intra Oral Camera
- LED Curing Light
- External Wet Suction
- High Speed Handpiece
- Low Speed Motor Handpiece
- Connection for External Suction
- Contra head & Straight Handpiece

## UPHOLSTERY COLOURS


Wide range of colours for stool & patient chair



VD# 01



VD# 02



VD# 03



VD# 04



VD# 05



VD# 06



VD# 07



VD# 08



VD# 09



VD# 10

Shown samples may appear different in structure and colour due to the printing process

future of healthcare

# VITANIX



+38664150661 | INFO@VITANIX.SI | WWW.VITANIX.SI

ADDRESS: Mlinska ulica 1 2000 Maribor Slovenia Europe

VITANIX D.o.o reserves the right to make any changes to the information, including specifications, without notice

