

DENTAL UNIT

DV PREMIERE-S



DV- PREMIERE - S

ACCURACY EFFICIENCY STABILITY
RELIABILITY & USER FRIENDLY



Touch pad for user Friendly operations

Stable design for your ultimate comfort

Integrated design to ensure its trustworthiness

Innovative control system to ensure precision
and effectiveness

Cuspidor Water Unit

Rotatable cuspidor spittoon & main chest
Automatic pressure clean water system

Assistant Element

The assistant operation rack is rotatable and is provided with an auxiliary control panel for controlling spittoon rinsing, mouthwash, and dental operating light, so as to facilitate the assistant operation
Equipped with safety protection for securer operation
Securer and easier operation for your beneficial experience
Assistant Element is easy to operate with trustworthy hands
Assistant Console have a holder left for easy upgrade, Dental Chair rinsing & warm water can be controlled by touch pad



OPERATION DENTAL LIGHT

Multi-axis rotation design to ensure that the light is irradiated at any Angle Longer life
Color temperature: 4000K yellow & 5000K white
LED Sensor light for no-shadow operation
Removable handle for easy disinfection
Stepless light adjustment
Illuminance: 6000-32000Lx

WORKING STOOLS

360° Swiveling backrest and seat
Manual hand release lever
Mounted on 5-castors



PATIENT DENTAL CHAIR

Right armrest is designed as a movable armrest for safely patient positioning

Further safety system under backrest and base plate ensures operation safety

User Friendly design for relaxing seating for comfort of patient in long treatments

The user can press any chair position button on the main control panel to stop any movement of the dental chair immediately

Safety switch is present when the backrest leans backward, it will stop movement immediately when encountering an obstacle

Chair safety system is present for dental chair and dental unit to work simultaneously

Synchronized movement of chair offer boundless comfort for patient

Height adjustment range of dental chair 415 mm-740 mm

Angle between backrest and horizontal plane: 0°-70°

Headrest adjustment range: 150 mm

High Quality dental Upholstery

Electric Chair Movement





Dentist Element

Air locking system for easy manipulation

Time and alarm Settings Priority function of first device Fault display

The instrument tray can bear a maximum load of 1.5 kg (uniform distribution)

Three different user interfaces can be provided the tray can be enlarged as required


Innovative & Latest instrument arm with two sets ball bearing ensures extraordinary performance

Rinsing time , Linkage function



DENTAL HANDPIECE

Optical fiber brightness can be adjusted at any level from 1 to 16



FOOT PEDAL

Lift and back motion of the chair
Free rotation of rotatable instruments
Operation light can be controlled via
Efficient foot pedal



Head Rest

Articulating Head Rest for fasy
Operations of Head Rest

STANDARDS ACCESSORIES

DENTIST ELEMENT

- 3-Way Syringe
- X-ray Film Viewer
- Small X-ray Film Viewer
- 3 outlets for turbine or air motor without Light

STOOL

- Gas Lift Cylinder
- Five castors
- Height adjustable lever

OPERATING LIGHT


- LED Light

WATER UNIT

- Moveable Cuspidor
- Water Bottle System
- Heat for cup Filler
- Venturi System

DENTAL CHAIR

- Trendelenburg Movement
- Multidimensional Control
- Safety protection Switches
- Seamless Upholstery
- Headrest 2 articulated
- Left & Right Armse
- Pre-set Positions



ASSISTANT ELEMENT

- Control Panel
- 3-Way Syringe
- Large Suction Outlet
- Small Suction Outlet
- Moveable Assistant Element

OPTIONS

- Air Compressor
- Intra Oral Camera
- Scaler
- LED Curing Light
- External Wet Suction
- High Speed Handpiece
- Low Speed Motor Handpiece
- Contra head & Straight Handpiece

UPHOLSTERY COLOURS


Wide range of colours for stool & patient chair



VD# 01



VD# 02



VD# 03



VD# 04



VD# 05



VD# 06



VD# 07



VD# 08



VD# 09



VD# 10

Shown samples may appear different in structure and colour due to the printing process

VITANIX

future
of
healthcare



+38664150661 | INFO@VITANIX.SI | WWW.VITANIX.SI

ADDRESS: Mlinska ulica 1 2000 Maribor Slovenia Europe

VITANIX D.o.o reserves the right to make any changes to the information, including specifications, without notice

