

AUTOCLAVE

VITA-STERL



VITA-STERL

ACCURACY EFFICIENCY STABILITY
RELIABILITY & USER FRIENDLY



Vita-Sterl is ultimate state of the art steam sterilizers
It is designed to achieve the highest level of sterilization
Microprocessor Controlled to perform complete sterilization
Monitoring system to check all the quality & programs
Temperature sensors ensure complete Monitoring of temperature



Accessibility

The 6.9" high-definition colour touch screen display, and its intuitive graphics, make the use of Vita-Sterl well-defined and trouble-free

Performance

High loading capacities and rapid sterilization times enable well-organized workflows

Safety

The integrated programs & quality interface provides supreme process reliability

Intelligent

Automatic identification of fault protection system

Technical Specification

The innovative Class B steam sterilizers for practice and clinic

- Tabletop Class B-Autoclave
- Sterilization 5 programs
- Rapid cycle programs
- Conductivity sensor
- Stainless-Steel chamber
- Temperature sensor
- Maximum Porous Load 2Kg
- Maximum Solid Load 5.5 Kg
- Built-in monitoring system
- Chamber Capacity: 23 Liter & 45 Liter
- Chamber Material Stainless Steel
- Tipple-pre and deep post vacuum
- Temperature Cycle 134°C in 40-44 min
- Temperature Cycle 121°C in 45-60 min
- Vacuum and Bowie & Dick test program
- Internal storage of last 45 cycles data
- Microprocessor Programmed & Smooth Sterilization
- Built-in Electric Jet Type Steam Generator S.S Material
- Automatic electrical door lock with safety door lock system
- 5 compatible trays
- Reservoirs for fresh and waste water are available for efficient sterilization
- Tray holder



Optional

- Distiller : Rapid water purification system
- Printer : Printing cycles
- Tray Holder : Tray holder for 3 or 5 shelves
- Tray : Trays for loading positioning

Model	Power	Voltage	Frequency	Tank Capacity	Appearance size	Inside Dimension	G · W
VITA- STERL 23	2000W	AC220V/110V ± 10%	50/60Hz ± 1Hz	23Litres	470X500X630	φ 249X450	59kg
VITA- STERL 45	3500W	AC220V/110V ± 10%	50/60Hz ± 1Hz	45Litres	540X580X800	φ 319X623	105kg

future of healthcare

VITANIX



+38664150661 | INFO@VITANIX.SI | WWW.VITANIX.SI

ADDRESS: Mlinska ulica 1 2000 Maribor Slovenia Europe

VITANIX D.o.o reserves the right to make any changes to the information, including specifications, without notice